

Anti-Cd74 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P04233-2

Item Number: CM0005260

Introduction

CamelBio supplies a rabbit monoclonal antibody targeting human CD74 (invariant chain, CLIP) validated for WB, IHC-P, IF/ICC, FC, and ELISA. Designed for MHC class II antigen presentation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD74 Rabbit mAb
Remarks/Alias	II; p33; CLIP; DHLAG; HLADG; Ia-GAMMA
Species	Human
Gene ID (Human)	972
Gene ID	972
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-232 of human CD74 (NP_004346.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P04233-2
Protein Weight	18kDa/24kDa/26kDa/31kDa/34kDa
Shipping	Ice bag