

Anti-Cd74 Rabbit Polyclonal Antibody For Wb, Elisa - P04441-2

Item Number: CM0023388

Introduction

Rabbit polyclonal antibody against mouse CD74, validated for WB and ELISA. Targets the invariant chain involved in MHC class II antigen presentation, expressed in thymus and lymph nodes.

[Learn More](#)

Parameter	Value
Product Name	CD74 Rabbit pAb
Remarks/Alias	Ii; CLIP; DHLAG; HLADG; Ia-GAMMA; CD74
Species	Mouse
GeneID	16149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-215 of mouse CD74 (NP_034675.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04441-2
Protein Weight	24kDa
Shipping	Ice bag