

Anti-Cd74 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P04233

Item Number: CM0025262



Introduction

High-quality anti-CD74 rabbit polyclonal antibody for WB, IF/ICC, IHC-P, ELISA. Detects human CD74 (P04233), cross-reacts with rat. Suitable for MHC class II chaperone, MIF receptor, and antiviral immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD74 Rabbit pAb
Remarks/Alias	II; p33; CLIP; DHLAG; HLADG; Ia-GAMMA; CD74
Species	Human
Gene ID (Human)	972
Gene ID	972
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-232 of human CD74 (NP_004346.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P04233
Protein Weight	34 kDa
Shipping	Ice bag