

## Anti-Cd74 Rabbit Polymab® - P04233-2 / P04441-2

Item Number: CM0012593

### Introduction

Recombinant rabbit monoclonal antibody targeting human and mouse CD74 (invariant chain). Validated for western blot, immunofluorescence, immunohistochemistry, and ELISA. Reacts with both p33 and p35 isoforms. P04233-2 / P04441-2.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD74 Rabbit PolymAb®                                                                                             |
| Remarks/Alias    | II; p33; CLIP; DHLAG; HLADG; Ia-GAMMA                                                                            |
| Species          | Human/Mouse                                                                                                      |
| Gene ID (Human)  | 972                                                                                                              |
| Gene ID (Mouse)  | 16149                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 73-232 of human CD74 (NP_004346.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, IF-P, IHC-P, ELISA                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | P04233-2/P04441-2                                                                                                |
| Protein Weight   | 24kDa/26kDa                                                                                                      |
| Shipping         | Ice bag                                                                                                          |