

Anti-Cd79A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11912

Item Number: CM0003404



Introduction

High-quality rabbit monoclonal antibody against human CD79a (MB-1, Ig-alpha) for WB, IHC-P, ELISA. Protein weight 25kDa. Research use for B-cell receptor signaling.

[Learn More](#)

Parameter	Value
Product Name	CD79a Rabbit mAb
Remarks/Alias	IGA; MB1; MB-1; IGA1pha; CD79a
Species	Human
GeneID (Human)	973
GeneID	973
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P11912
Protein Weight	25kDa
Shipping	Ice bag