

Anti-Cd79A Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11912

Item Number: CM0015145



Introduction

Rabbit polyclonal antibody raised against recombinant CD79a (aa 165-225). Validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Detects ~25kDa protein. Ideal for B-cell receptor signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD79a Rabbit pAb
Remarks/Alias	IGA; MB1; MB-1; IGA1pha; CD79a
Species	Human
GeneID (Human)	973
GeneID	973
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 165-225 of human CD79a (NP_001774.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11912
Protein Weight	25kDa
Shipping	Ice bag