

Anti-Cd79B Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P40259

Item Number: CM0006369



Introduction

Recombinant anti-CD79B rabbit monoclonal antibody validated for WB, IF/ICC, FC, ELISA. Detects human CD79B (14/26 kDa) derived from recombinant immunogen. Suitable for B-cell receptor signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD79B Rabbit mAb
Remarks/Alias	CD79B; AGM6; B29; IGB; CD79b molecule
Species	Human
GeneID (Human)	974
GeneID	974
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79B (NP_000617.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P40259
Protein Weight	14kDa/26kDa
Shipping	Ice bag