

Anti-Cd79B Rabbit Monoclonal Antibody (B29, Ig-Beta) For Wb, If, Fc, Elisa - P15530

Item Number: CM0006866



Introduction

Rabbit monoclonal antibody specific to mouse CD79b (B29, Ig-beta), confirmed for WB, IF, FC, ELISA. Detects ~25 kDa band; cross-reacts with rat. Ideal for B-cell receptor studies.

[Learn More](#)

Parameter	Value
Product Name	CD79b Rabbit mAb
Remarks/Alias	B29; Igb; Igbeta; Ig-beta
Species	Mouse
GeneID (Human)	974
GeneID	15985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse CD79b (NP_032365.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, FC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P15530
Protein Weight	25kDa
Shipping	Ice bag