

Anti-Cd79B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P40259

Item Number: CM0018904



Introduction

High-quality rabbit polyclonal antibody targeting human CD79B (B29, Igbeta, AGM6). Validated for WB, IHC-P, ELISA; cross-reacts with mouse, rat. Ideal for B-cell receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD79B Rabbit pAb
Remarks/Alias	B29; IGB; AGM6; Igbeta; CD79B
Species	Human
GeneID (Human)	974
GeneID	974
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 130-229 of human CD79B (NP_000617.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40259
Protein Weight	26kDa
Shipping	Ice bag