

Anti-Cd79B Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P40259

Item Number: CM0006428



Introduction

CD79B Rabbit PolymAb® is a monoclonal antibody targeting human CD79B, validated in WB, IF-P, IHC-P, ELISA, and cross-reacts with mouse and rat. Ideal for B-cell receptor signaling and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD79B Rabbit PolymAb®
Remarks/Alias	B29; IGB; AGM6; Igbeta
Species	Human
GeneID (Human)	974
GeneID	974
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 130-229 of human CD79B (NP_000617.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40259
Protein Weight	26kDa
Shipping	Ice bag