

Anti-Cd80 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P33681

Item Number: CM0006939



Introduction

High-quality Anti-CD80 Rabbit Monoclonal Antibody validated for Western Blot, Flow Cytometry, and ELISA. Specific to human CD80, a critical T-cell costimulatory molecule. Suitable for immunology and immuno-oncology research.

[Learn More](#)

Parameter	Specification
Product Name	CD80 Rabbit mAb
Remarks/Alias	CD80; B7; B7-1; B7.1; BB1; CD28LG; CD28LG1; LAB7; CD80 molecule
Species	Human
GeneID	941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-242 of human CD80 (NP_005182.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P33681
Protein Weight	18kDa/29kDa/33kDa
Shipping	Ice bag