

## Anti-Cd80 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q00609

Item Number: CM0007285



### Introduction

Rabbit monoclonal antibody against mouse CD80, validated for WB, FC, ELISA. Recognizes human and mouse CD80. Suitable for T-cell activation and immune checkpoint research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD80 Rabbit mAb                                                                                                  |
| Remarks/Alias    | B71; Ly53; TSA1; Cd28I; Ly-53; MIC17; CD80                                                                       |
| Species          | Mouse                                                                                                            |
| GeneID (Human)   | 941                                                                                                              |
| GeneID           | 12519                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 38-245 of mouse CD80 (NP_033985.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | WB, FC, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | Q00609                                                                                                           |
| Protein Weight   | 34 kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |