

Anti-Cd80 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P33681

Item Number: CM0024289



Introduction

High-quality anti-CD80 rabbit polyclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Recognizes human, mouse, and rat CD80 (B7-1), a key costimulatory molecule for T-cell activation. Ideal for immunology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD80 Rabbit pAb
Remarks/Alias	B7, BB1, B7-1, B7.1, LAB7, CD28LG, CD28LG1, CD80
Species	Human
GeneID (Human)	941
GeneID	941
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD80 (NP_005182.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33681
Protein Weight	33kDa
Shipping	Ice bag