


Anti-Cd81 Rabbit Polyclonal Antibody For Wb, Elisa - P35762

Item Number: CM0020856



Introduction

High-quality rabbit polyclonal antibody against mouse CD81 (TAPA-1, Tspan28), validated for Western blot and ELISA, cross-reacts with rat. Suitable for immune cell signaling, tetraspanin biology, and host-pathogen interaction studies.

[Learn More](#)

Parameter	Value
Product Name	CD81 Rabbit pAb
Remarks/Alias	Tapa1; Tapa-1; Tspan28; CD81
Species	Mouse
Gene ID (Human)	975
Gene ID (Mouse)	12520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 106-211 of mouse CD81 (NP_598416.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35762
Protein Weight	25 kDa
Shipping	Ice bag