


Anti-Cd81 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P60033

Item Number: CM0025985



Introduction

Anti-CD81 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. CD81 is a tetraspanin involved in B cell signaling, HCV entry, and immune synapse formation. Ideal for immunology, infectious disease, and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CD81 Rabbit pAb
Remarks/Alias	S5.7; CVID6; TAPA1; TSPAN28; CD81
Species	Human
Gene ID (Human)	975
Gene ID	975
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 114-199 of human CD81 (NP_004347.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60033
Protein Weight	26kDa
Shipping	Ice bag