


Anti-Cd81 Rabbit Polymab® For Wb, Ihc-P, Elisa - P60033

Item Number: CM0004298



Introduction

Rabbit recombinant monoclonal antibody to human CD81 (TAPA1/TSPAN28). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying tetraspanin-enriched microdomains, B cell signaling, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	CD81 Rabbit PolymAb®
Remarks/Alias	S5.7; CVID6; TAPA1; TSPAN28
Species	Human
Gene ID (Human)	975
Gene ID	975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 114-199 of human CD81 (NP_004347.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60033
Protein Weight	26kDa
Shipping	Ice bag