

Anti-Cd82 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27701

Item Number: CM0020131



Introduction

High-quality rabbit polyclonal antibody targeting CD82 (KAI1/TSPAN27), a metastasis suppressor and tetraspanin. Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for cancer and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	CD82 Rabbit pAb
Remarks/Alias	R2, 4F9, C33, IA4, ST6, GR15, KAI1, SAR2, TSPAN27, CD82
Species	Human
GeneID (Human)	3732
GeneID	3732
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human CD82 (NP_002222.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27701
Protein Weight	30kDa
Shipping	Ice bag