

Anti-Cd82 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P27701

Item Number: CM0007884



Introduction

CD82 rabbit monoclonal antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Detects human and mouse CD82 metastasis suppressor, a tetraspanin involved in cell adhesion, migration, and signaling. Ideal for cancer and immunology research.

[Learn More](#)

Parameter	Description
Product Name	CD82 Rabbit mAb
Remarks/Alias	R2, 4F9, C33, IA4, ST6, GR15, KAI1, SAR2, TSPAN27, CD82
Species	Human
Gene ID	3732
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 168-267 of human CD82 (P27701)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27701
Protein Weight	30kDa
Shipping	Ice bag