

Anti-Cd83 Rabbit Monoclonal Antibody For Wb, Elisa - Q01151

Item Number: CM0007677



Introduction

Rabbit monoclonal antibody against human CD83, a transmembrane glycoprotein on immune cells, suitable for WB and ELISA research applications.

[Learn More](#)

Parameter	Value
Product Name	CD83 Rabbit mAb
Remarks/Alias	BL11; HB15; CD83
Species	Human
GeneID (human)	9308
GeneID	9308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-143 of human CD83 (NP_004224.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q01151
Protein Weight	23 kDa
Shipping	Ice bag