

Anti-Cd83 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q01151

Item Number: CM0007831



Introduction

Rabbit monoclonal antibody against human CD83, validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for immune response and dendritic cell research.

[Learn More](#)

Parameter	Value
Product Name	CD83 Rabbit mAb
Remarks/Alias	BL11; HB15; CD83
Species	Human
Gene ID (Human)	9308
Gene ID	9308
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD83 (Q01151).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01151
Protein Weight	23kDa
Shipping	Ice bag