


Anti-Cd84 Monoclonal Antibody For Wb And Elisa - Q9Uib8

Item Number: CM0007321



Introduction

Rabbit monoclonal antibody against human CD84 (SLAMF5), validated for WB and ELISA. Cross-reacts with mouse and rat. Targets extracellular domain; suitable for immune cell signaling and functional studies.

[Learn More](#)

Parameter	Value
Product Name	CD84 Rabbit mAb
Remarks/Alias	LY9B; hCD84; mCD84; SLAMF5; CD84
Species	Human
Gene ID	8832
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 246-345 of human CD84 (Q9UIB8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIB8
Protein Weight	39kDa
Shipping	Ice bag