


Anti-Cd84 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q9Uib8

Item Number: CM0004272



Introduction

Rabbit monoclonal antibody targeting CD84 (SLAMF5), designed for WB, IF/ICC, FC, and ELISA applications. Reacts specifically with human CD84, a 39 kDa immune cell receptor important for signaling and hematopoietic studies.

[Learn More](#)

Parameter	Value
Product Name	CD84 Rabbit mAb
Remarks/Alias	LY9B; hCD84; mCD84; SLAMF5; CD84
Species	Human
Gene ID	8832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-225 of human CD84 (NP_001171808.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	Q9UIB8
Protein Weight	39kDa
Shipping	Ice bag