


Anti-Cd84 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uib8

Item Number: CM0026690



Introduction

High-quality rabbit polyclonal antibody targeting human CD84 (SLAMF5). Validated for Western blot, immunofluorescence, and ELISA. Ideal for studying immune cell activation, SLAM family signaling, and hematopoietic progenitor cell biology.

[Learn More](#)

Parameter	Value
Product Name	CD84 Rabbit pAb
Remarks/Alias	LY9B; hCD84; mCD84; SLAMF5; CD84
Species	Human
GeneID (Human)	8832
GeneID	8832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-225 of human CD84 (NP_003865.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UIB8
Protein Weight	39kDa
Shipping	Ice bag