


Anti-Cd85j Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nhl6

Item Number: CM0024525



Introduction

High-quality rabbit polyclonal antibody targeting CD85j (LILRB1), validated for Western blot, IF/ICC, and ELISA. Reacts with human samples. An essential tool for immune regulation and inhibitory receptor research.

[Learn More](#)

Parameter	Value
Product Name	CD85j Rabbit pAb
Remarks/Alias	ILT2; LIR1; MIR7; PIRB; CD85j; ILT-2; LIR-1; MIR-7; PIR-B; CD85j
Species	Human
Gene ID (Human)	10859
Gene ID	10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-458 of human CD85j (NP_006660.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8NHL6
Protein Weight	71kDa
Shipping	Ice bag