


Anti-Cd86 Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P42081

Item Number: CM0002784



Introduction

This rabbit monoclonal antibody targets human CD86 (B7-2, LAB72), a costimulatory molecule on antigen-presenting cells. Validated for WB, IF/ICC, FC, and ELISA, it supports immune activation and checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	CD86 Rabbit mAb
Remarks/Alias	B70; B7-2; B7.2; LAB72; CD28LG2; CD86
Species	Human
Gene ID (Human)	942
Gene ID	942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-241 of human CD86 (NP_787058.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P42081
Protein Weight	38kDa
Shipping	Ice bag