


Anti-Cd86 Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, If-P, Elisa - P42081

Item Number: CM0023954



Introduction

CD86 rabbit polyclonal antibody validated for WB, IP, IF/ICC, IF-P, ELISA. Ideal for studying T-cell costimulation, CD28/CTLA4 pathways, and viral entry mechanisms. Reacts with human, mouse, rat CD86.

[Learn More](#)

Parameter	Value
Product Name	CD86 Rabbit pAb
Remarks/Alias	B70; B7-2; B7.2; LAB72; CD28LG2; CD86
Species	Human
GeneID (Human)	942
GeneID	942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human CD86 (NP_787058.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42081
Protein Weight	38kDa
Shipping	Ice bag