


Anti-Cd89 Monoclonal Antibody For Flow Cytometry - P24071

Item Number: CM0009013



Introduction

-conjugated rabbit anti-human CD89 monoclonal antibody, validated for flow cytometry with reactivity against human samples. Target: immunoglobulin alpha Fc receptor (FCAR, FcalphaRI, 32 kDa).

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD89 mAb
Remarks/Alias	CD89; FcalphaR; FcalphaRI; CTB-61M7.2
Species	Human
Gene ID (Human)	2204
Gene ID	2204
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24071
Protein Weight	32kDa
Shipping	Ice bag