


Anti-Cd89 Rabbit Monoclonal Antibody For Flow Cytometry - P24071

Item Number: CM0009774



Introduction

Rabbit monoclonal antibody against human CD89 (FcalphaR), validated for flow cytometry. Reacts with human samples, molecular weight 32kDa. Ideal for studying IgA-mediated immune responses. Shipped on ice bags.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD89 mAb
Remarks/Alias	CD89; FcalphaR; FcalphaRI; CTB-61M7.2
Species	Human
Gene ID	2204
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-227 of human CD89 (NP_001991.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24071
Protein Weight	32kDa
Shipping	Ice bag