

Anti-Cd8A Monoclonal Antibody For Fc - P07725

Item Number: CM0008010

Introduction

ABflo 647-conjugated rabbit monoclonal antibody targeting rat CD8a for flow cytometry (FC). Recombinant immunogen, reactive with rat. Ideal for CD8+ T-cell analysis, immune cell phenotyping, and cytotoxic lymphocyte research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Rat CD8a mAb
Remarks/Alias	T8; Lyt2; Ly-2; CD8α
Species	Rat
GeneID (Human)	925
GeneID	24930
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-189 of rat CD8a (NP_113726.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	P07725
Protein Weight	32kd
Shipping	Ice bag