


Anti-Cd8A Polyclonal Antibody For Wb, Ihc-P, Elisa - P01732

Item Number: CM0018611



Introduction

Robust rabbit polyclonal antibody against human CD8A, validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse CD8A. Suitable for T-cell and immune response research.

[Learn More](#)

Parameter	Specification
Product Name	CD8A Rabbit pAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha; CD8A
Species	Human
GeneID (Human)	925
GeneID	925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01732
Protein Weight	26kDa
Shipping	Ice bag