


Anti-Cd8A Rabbit Polyclonal Antibody For Wb And Elisa - P01731

Item Number: CM0018925



Introduction

Rabbit polyclonal antibody against mouse CD8A (Lyt-2), validated for Western blot and ELISA. Recognizes the 27 kDa T-cell coreceptor, ideal for immune response and cytotoxic T-cell research.

[Learn More](#)

Parameter	Value
Product Name	CD8A Rabbit pAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2; CD8A
Species	Mouse
GeneID (Human)	925
GeneID (Mouse)	12525
Immunogen	Recombinant protein of mouse CD8A
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27 kDa
Shipping	Ice bag