


Anti-Cd8A Rabbit Polymab® For Wb, Ihc-P, Mihc - P01732 / P01731

Item Number: CM0004007



Introduction

Recombinant rabbit monoclonal antibody targeting human/mouse CD8A. Validated for WB, IF-P, IHC-P, mIHC, and ELISA. Detects CD8A across human, mouse, rat, and cynomolgus monkey. Essential for T-cell immunology and immuno-oncology studies.

[Learn More](#)

Parameter	Details
Product Name	CD8A Rabbit PolymAb®
Remarks/Alias	CD8A;CD8;Leu2;MAL;p32
Species	Human/Mouse
Gene ID (Human)	925
Gene ID	925/12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat, Cynomolgus monkey
SWISS	P01732/P01731
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag