

Anti-Cd8B Polyclonal Antibody For Wb, Elisa - P10966

Item Number: CM0029451

Introduction

High-quality rabbit polyclonal antibody against human CD8B (CD8beta). Validated for Western blot and ELISA, cross-reactive with human, mouse, and rat. Recognizes the T-cell surface glycoprotein CD8 beta chain, essential for MHC class I interaction and CTL activation.

[Learn More](#)

Parameter	Value
Product Name	CD8B Rabbit pAb
Remarks/Alias	LY3; P37; LEU2; LYT3; Ly-3; CD8B1; CD8beta; CD8B
Species	Human
GeneID (Human)	926
GeneID	926
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 22-170 of human CD8B (NP_757362.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10966
Protein Weight	24kDa
Shipping	Ice bag