

Anti-Cd8B Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P10966

Item Number: CM0006814

Introduction

Rabbit monoclonal antibody against human CD8 β (CD8B, LYT3) validated for Western blotting, flow cytometry, and ELISA. Suitable for T-cell immunology and immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	CD8 β Rabbit mAb
Remarks/Alias	LY3, P37, LEU2, LYT3, Ly-3, CD8B1, CD8beta, CD8 β
Species	Human
GenelD (human)	926
GenelD	926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-170 of human CD8B (NP_004922.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P10966
Protein Weight	24kDa
Shipping	Ice bag