


Anti-Cd9 Rabbit Monoclonal Antibody For Wb, If, Ihc-P, Fc, Elisa - P21926

Item Number: CM0007016



Introduction

CD9 rabbit monoclonal antibody validated for Western blot, immunofluorescence, immunohistochemistry, flow cytometry, and ELISA; detects human CD9 tetraspanin involved in cell adhesion, sperm-egg fusion, and platelet activation. For research use.

[Learn More](#)

Parameter	Value
Product Name	CD9 Rabbit mAb
Remarks/Alias	MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29; CD9
Species	Human
Gene ID (Human)	928
Gene ID	928
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P21926
Protein Weight	25 kDa
Shipping	Ice bag