


Anti-Cd9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P21926

Item Number: CM0019823



Introduction

High-quality CD9 Rabbit pAb (P21926) for WB, IHC-P, ELISA. Validated in human, mouse, rat samples. Sourced from rabbit, immunogen a.a. 112-195. Reliable for cell adhesion, platelet activation, and sperm-egg fusion research.

[Learn More](#)

Parameter	Value
Product Name	CD9 Rabbit pAb
Remarks/Alias	MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29; CD9
Species	Human
GeneID (Human)	928
GeneID	928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-195 of human CD9 (NP_001760.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21926
Protein Weight	25kDa
Shipping	Ice bag