

Anti-Cd90.1 Rabbit Polyclonal Antibody For Wb, Elisa - P01831

Item Number: CM0020585

Introduction

Rabbit polyclonal antibody against mouse CD90.1 (Thy1.1) for WB and ELISA applications. Targets GPI-anchored Thy-1 membrane glycoprotein involved in cell adhesion and synaptogenesis.

[Learn More](#)

Parameter	Value
Product Name	CD90.1 Rabbit pAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2; CD90.1
Species	Mouse
Gene ID (Human)	7070
Gene ID (Mouse)	21838
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-120 of mouse CD90.1 (NP_033408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag