

Anti-Cd90.2 Rabbit Monoclonal Antibody For Elisa, Fc - P01831

Item Number: CM0011743

Introduction

CD90.2 Rabbit monoclonal antibody (mAb) validated for ELISA and flow cytometry (FC). Recognizes mouse Thy-1 membrane glycoprotein (P01831). Ideal for immunology and neurobiology research applications.

[Learn More](#)

Parameter	Value
Product Name	CD90.2 Rabbit mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
GeneID (Human)	7070
GeneID (Mouse)	21838
Immunogen	Recombinant protein of mouse CD90.2.
Source	Rabbit
Category	Monoclonal Antibodies
Application	ELISA, FC
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag