

Anti-Cd90/Thy1 Monoclonal Antibody For Fc - P04216

Item Number: CM0004225

Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD90/Thy1, validated for flow cytometry (FC). Ideal for cell surface CD90 expression analysis in immunology and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD90/Thy1 mAb
Remarks/Alias	CD90; CDw90
Species	Human
GeneID (Human)	7070
GeneID	7070
Immunogen	Recombinant protein (recombinant fusion protein containing amino acids 20-130 of human CD90/Thy1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag