

Anti-Cd90/Thy1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P04216

Item Number: CM0002795

Introduction

Rabbit monoclonal antibody targeting CD90/Thy1 (UniProt P04216) for western blot, immunohistochemistry, and ELISA. Detects human and rat Thy-1 membrane glycoprotein involved in cell-cell interactions and synaptogenesis.

[Learn More](#)

Parameter	Value
Product Name	CD90/Thy1 Rabbit mAb
Remarks/Alias	CD90; CDw90; CD90/Thy1
Species	Human
Gene ID (Human)	7070
Gene ID	7070
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD90/Thy1 (P04216)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag