

Anti-Cd90/Thy1 Monoclonal Antibody (Kd Validated) For Wb, Fc, Elisa - P04216

Item Number: CM0007282

Introduction

High-quality recombinant rabbit monoclonal antibody against human CD90/Thy1 (P04216), KD validated. Ideal for western blot, flow cytometry, and ELISA. Investigates cell-cell interactions and synaptogenesis.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CD90/Thy1 Rabbit mAb
Remarks/Alias	CD90; CDw90; [KO/KD Validated] CD90/Thy1
Species	Human
GeneID (Human)	7070
GeneID	7070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CD90/Thy1 (NP_006279.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag