

Anti-Cd90/Thy1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01831

Item Number: CM0008305

Introduction

High-quality CD90/Thy1 rabbit monoclonal antibody for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Recognizes human, mouse, and rat Thy-1 antigen. Suitable for neuroscience and cell membrane studies.

[Learn More](#)

Parameter	Value
Product Name	CD90/Thy1 Rabbit mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
GeneID (Human)	7070
GeneID	21838
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of mouse CD90/Thy1 (NP_033408.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag