

Anti-Cd90/Thy1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P04216

Item Number: CM0021221

Introduction

Rabbit polyclonal antibody against human CD90/Thy1, validated for WB, IF/ICC, and ELISA. Detects endogenous Thy-1 membrane glycoprotein in human, mouse, and rat samples. Ideal for neuroscience and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD90/Thy1 Rabbit pAb
Remarks/Alias	CD90; CDw90; CD90/Thy1
Species	Human
GeneID (human)	7070
GeneID	7070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CD90/Thy1 (NP_006279.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag