

Anti-Cd93 Rabbit Monoclonal Antibody For If/Icc, Fc, Elisa - O89103

Item Number: CM0004083

Introduction

Rabbit monoclonal antibody targeting mouse CD93 (C1q receptor, Ly68) for IF/ICC, flow cytometry, and ELISA. Cross-reacts with human CD93. Ideal for studying phagocytosis, inflammation, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	CD93 Rabbit mAb
Remarks/Alias	Ly68; AA4.1; C1qr1; C1qrp; 6030404G09Rik
Species	Mouse
GeneID (Human)	22918
GeneID	17064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-181 of mouse CD93 (NP_034870.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O89103
Protein Weight	69kDa
Shipping	Ice bag