

## Anti-Cd94 Monoclonal Antibody For Fc - Q13241

Item Number: CM0005899



### Introduction

Rabbit monoclonal antibody targeting human CD94 (KLRD1), conjugated with . Validated for flow cytometry, suitable for NK cell and T cell immune monitoring research.

[Learn More](#)

| Parameter        | Value                                          |
|------------------|------------------------------------------------|
| Product Name     | Rabbit anti-Human CD94 mAb                     |
| Remarks/Alias    | KLRD1; CD94; natural killer cells antigen CD94 |
| Species          | Human                                          |
| GeneID           | 3824                                           |
| Immunogen        | Recombinant protein of human CD94              |
| Source           | Rabbit                                         |
| Category         | Monoclonal Antibodies                          |
| Application      | FC                                             |
| Cross Reactivity | Human                                          |
| SWISS            | Q13241                                         |
| Protein Weight   | 17kDa/20kDa                                    |
| Shipping         | Ice bag                                        |