


Anti-Cd95 Rabbit Monoclonal Antibody For Elisa, Fc - P25446

Item Number: CM0006041



Introduction

Recombinant rabbit monoclonal antibody targeting mouse CD95/FAS (aa 22-169). Validated for ELISA and flow cytometry. Recognizes the apoptosis-mediating receptor FAS in mouse tissues. Ideal for apoptosis and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD95 Rabbit mAb
Remarks/Alias	lpr; APO1; APT1; CD95; TNFR6; Tnfrsf6
Species	Mouse
GeneID (human)	355
GeneID	14102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-169 of CD95 (NP_032013.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	ELISA, FC
Cross Reactivity	Mouse
SWISS	P25446
Protein Weight	37kDa
Shipping	Ice bag