


Anti-Cd95/Fas Monoclonal Antibody () For Flow Cytometry - P25446

Item Number: CM0008887



Introduction

Rabbit monoclonal antibody against mouse CD95/FAS, directly conjugated with for flow cytometry. Suitable for apoptosis research in thymus, liver, lung, and heart tissues.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD95/FAS mAb
Remarks/Alias	lpr; APO1; APT1; CD95; TNFR6; Tnfrsf6
Species	Mouse
Gene ID (Human)	355
Gene ID (Mouse)	14102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-169 of mouse CD95/FAS. (NP_032013.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25446
Protein Weight	37KDa
Shipping	Ice bag