

Anti-Cd95/Fas Monoclonal Antibody For Wb, If/Icc, Elisa - P25445

Item Number: CM0005642



Introduction

Recombinant rabbit monoclonal antibody targeting human CD95/FAS (UniProt P25445). Validated for Western blot, IF/ICC, and ELISA, supporting apoptosis and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD95/FAS Rabbit mAb
Remarks/Alias	APT1; CD95; FAS1; APO-1; FASTM; ALPS1A; TNFRSF6
Species	Human
GeneID (Human)	355
GeneID	355
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 26-173 of human CD95/FAS (NP_000034.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P25445
Protein Weight	9-38 kDa
Shipping	Ice bag