

Anti-Cd96 Polyclonal Antibody For Wb, Elisa - P40200

Item Number: CM0024294



Introduction

High-quality rabbit polyclonal antibody against human CD96 (TACTILE), validated for WB and ELISA. Detects CD96 in human, mouse, and rat samples. Ideal for studying T cell and NK cell adhesion and late-phase immune responses.

[Learn More](#)

Parameter	Value
Product Name	CD96 Rabbit pAb
Remarks/Alias	TACTILE; CD96
Species	Human
GeneID (Human)	10225
GeneID	10225
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-370 of human CD96 (NP_937839.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40200
Protein Weight	66kDa
Shipping	Ice bag