

Anti-Cd96 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P40200

Item Number: CM0007589



Introduction

CD96 rabbit monoclonal antibody validated for Western blot, flow cytometry, and ELISA. Recognizes human CD96, a T and NK cell adhesion receptor involved in late immune responses. Suitable for research on immune cell interactions.

[Learn More](#)

Parameter	Value
Product Name	CD96 Rabbit mAb
Remarks/Alias	Cell surface antigen CD96; CD96
Species	Human
GeneID (Human)	10225
GeneID	10225
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-519 of human CD96 (NP_937839.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P40200
Protein Weight	65kDa
Shipping	Ice bag